PORTUGUÊS

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INSTRUMENTAL CIRÚRGICO CORTANTE

Definição:

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ENGLISH

DEFINITION

INVASIVE DEVICES Rule / Classification 6/1 (according to EC Directive 93/42 / ECC)

PRESENTATION

This device is designed to specialized procedure, which should be performed by Dental Implant licensed professionals. For better results, use the product with the knowledge of the proper techniques. Always use them in appropriate conditions, including surgical room.

DESCRIPTION

The Cutting Instrumental is manufactured in surgical stainless steel. In its active end there are cylindrical profile and cutting geometry compatible with the DSP implants. The Twist Drill receives the application of a dark color thin carbon film on the cutting surface. This film has as main purposes:

- reduction of the heat generated during osteotomy
- friction reduction between the Drill and the bone
- resistance increase to wear • resistance increase to oxidation.

APPLICATIONS

This product is a surgical instrument used for drilling the bone tissue, with the purpose of implant bed preparation for insertion of the implant chosen.

DRILL LAUNCHES: Consisting of extremely sharp long position indicated to start the cortical bone drilling. Its use is not intended to drill in depth according to the size of the desired implant, but rather to break the cortical bone strength, being sufficient its deepening until the

cancellous bone is reached. It is the initial drill to all perforations, and in the case of tooth extractions with immediate implant placement can not be used.

DRILL SCREW: Suitable for drilling the bone surface, allowing the installation of cylindrical bodies implants. Must cortical bone must first be broken using the Drill lance, as directed above.

CONE DRILL: It has the same function as twist drill, but is recommended for bone drilling for installation of tapered bodies implants.

DRILL PILOT: It's an intermediate drill. It has a blunt tip, corresponding to the width immediately below drill and a cutting portion just above the corresponding thickness the following drill drilling. This feature allows the dentist to perform an initial drilling more cortical portion of the bone, without deviating from the initially planned drilling.

DRILL COUNTERSINK: Suitable for installation of implants with a diameter of platform above the implant body (implants also called expanded platform).

DRILL trephine: It must be used under the same conditions of drill bits. Intended to two situations: 1) obtaining graft circular section blocks, for example, the processing or maxillary tuberosity; 2) removal of misplaced implants or who have suffered damage to the internal thread. It has very cutting surface and should be used very carefully in areas near vital structures.

LIFT BONE: It has a sharp concave long position to charge and push the superior maxillary sinus floor bone. Use for socket preparation for the installation of implants in the posterior maxilla with low bone height.

BONE EXPANDER / RATCHET: Used when the thickness of the bone in the maxilla is inadequate or poor quality requiring compression and prior bone expansion to the implant platform. For the use of 3.8 mm measured expander, a milling cutter employing 2.0 to carry out the initial perforation, then followed by insertion of the 3.8 mm drilling expander to expand and compress the bone plate this particular measure. For the use of 4.3 mm measured expander, a milling cutter employing 2.0 to carry out the initial perforation, then a milling and below 2.15 or 2.5 to 4.3 mm expander insert in the bore to expand and compress the bone plate this specific measure. For use as 5.0 mm expander, a milling cutter employing 2.0 to carry out the initial perforation, then 2.5 cutter, then the cutter and follows 3.0 to 5.0 mm expander insert in the bore to expand and compress the table bone that particular measure. Obs.: Bone very thin, consisting of cortical, without the presence of cancellous bone, does not respond adequately to the expansion technique. It is recommended alternative technique of bone grafting.

MALE THREAD: Used for making bone contour for introducing Threaded implant when used in bones type I and II.

SCALPEL GINGIVAL: circular scalpel to cut gingival for inserting an implant or opening in second surgical procedure.

CONTRAINDICATIONS

This product does not present contraindications, once they are properly used for the indicated purposes.

HANDLING

Fit the Drill in the Contra-Angle. Place it in the specific rotation, according to instructions for use of the implant chosen. Perform the drilling with complete insertion and removal movements of the Drill. During drilling, the pressure should not be excessive and, an abundant irrigation, manual or combined, using the handpiece pump is necessary. Note 1: Do not interrupt the handpiece rotation with the Drill within the surgical cavity, because this can make difficult the removal or cause the Drill break.

SANITATION

The surgical instruments should be correctly sanitized after each use. Do not leave the product in contact with humidity more than the time necessary for the sanitation.

Proceed as follows:

- 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer).
- 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes.
- 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon

brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air.

5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent – 1st step – and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6th step: If applicable, place the product in an appropriate case duly sanitized and dry.

7th step: Select the packing in accordance to the sterilization process. Use preferably surgical grade paper self-sealing packing with laminated film.

FORMS OF PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 – Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave.

Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

PRECAUTIONS

- Do not use the product if the packing is damaged.
- This product should be sterile.
- The non-replacement of the Drills, as recommended by the manufacturer, can generate improper bone heating, compromising the osteointegration success.

• Do not regard.

• Do not use solutions with ammonia, acids or hydrogen peroxide, because they can damage the Drill.

• Do not use antirust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that will have cleaning solution wastes, because these procedures cause oxidation.

• The use of steel brushes is inadvisable. We recommend the use of nylon brushes.

• The improper planning can jeopardize the performance of the implant/prosthesis set causing failures to the system, such as loss or break of the Implant, loosening or fracture of the prosthetic screws.

• Better results are obtained with the use of the DSP products sequence.

• The use of instruments and/or prosthetic abutments of other manufacturers do not guarantee the perfect function of the DSP and exempts any guarantee of the product.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

POSTOPERATIVE CARE AND MAINTENANCE

Guide the patient as for the need to perform a professional monitoring after the surgery and to obey the guidelines related to care, hygiene and prescription of drugs. These guidelines are professional's responsibility.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIALS

All consumption products and materials used in the surgery for the insertion of dental implants can endanger the health of those who handle, after the use. Before disposing them into the environment, we recommend to refer to and comply with the effective legislation.

SHELF-LIFE: Suitable on the label.

LIFE CYCLE

This product has a life cycle of up to 20 drillings for bone qualities I and II, and up to 30 drillings for bone qualities III and IV, since its use conditions are respected, including adequate irrigation and rotation. Apart of the life cycle period mentioned, the professional should check the sharpening conditions of the Drills. Drills with improper cutting conditions should be disposed.

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Proteger contra luz
Keep away from sun

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