# DSP *white*

CLINIC 35% CALCIUM

EN | ENGLISH

PT | PORTUGUÊS

ES | ESPAÑOL

DSP  
*white*

# DSP white

CLINIC 35% CALCIUM

ENGLISH

To access the IFU go to  
dspbiomedical.com/ifu  
with code: 0.5005

TRADE NAME: DSP White Clinic 35% - Calcium  
TECHNICAL NAME: Dental Whitening Agents

Please read all information of this instruction for use carefully before applying the product.

DENTAL WHITENER FOR CLINICAL USE.  
FOR HOME APPLICATION UNDER THE SUPERVISION OF A DENTAL SURGEON  
SALE UNDER PRESCRIPTION (ACCORDING TO RDC NO. 6 OF 2015)

Keep for consultation until the total consumption of the product.  
Keep this instruction until there is no more interaction of the product with your patient.

## CONTENT:

- KIT**
- 6 Thickener syringes,
  - 6 Hydrogen Peroxide syringes,
  - 6 syringe connectors,
  - 1 bottle of Peroxide Neutralizer,
  - 1 DSP Dam syringe with 2g (with applicator tips),
  - 1 instruction manual for the professional.

## MINI KIT

- 1 Thickener syringes,
- 1 Hydrogen Peroxide syringes,
- 1 syringe connector,
- 1 DSP Dam syringe 1g (with applicator tip),
- 1 instruction manual for the professional.

## PRODUCT DESCRIPTION:

DSP White Clinic 35% - Calcium is a professional dental whitening system based on Hydrogen Peroxide 35% for whitening of vital teeth, for exclusive use in the dental clinic bringing the highest technology and new trends in world dentistry in terms of conservative techniques and applications and with more security for whitening in the clinic. And it should always be used with an efficient isolation system for gingival protection (DSP Dam light-cured gingival barrier). After the homogeneous mixing of the two phases, they form a 35% Hydrogen Peroxide gel, with controlled pH and containing desensitizing agents in its formula that maximize the patient's comfort.

The action of whitening agents is related to their penetration into dental structures (enamel and dentin) reacting with the pigmented molecules, oxidizing them and making them simpler and less pigmented, causing the pigments incorporated into the tissues to be totally or partially eliminated from the dental structure.

The auto mixing system provided by the syringes ensures the ideal concentration of peroxide and its effectiveness with each application. Once mixed, phases A and B are ready for use.

## PRODUCT INDICATION:

DSP White Clinic 35% - Calcium is indicated for vital teeth with extrinsic color changes of a chronic nature, not susceptible to removal by conventional procedures (prophylaxis) and for teeth with less severe intrinsic color changes (physiological darkening by dentin deposition) traumas, etc. In the most severe intrinsic color changes such as fluorosis, severe tetracycline staining, although effective, the results of whitening are less predictable. The use of whiteners is contraindicated in cases of imperfect amelogenesis and dentinogenesis.

DSP White Clinic 35% - Calcium can also be indicated as a complementary technique for whitening non-vital teeth associated with vitalized teeth whitening.

## COMPOSITION:

Hydrogen Peroxide 35%, Thickeners, neutralizers, dye mixture, Silicas, Calcium, glycols and Deionized Water.

## INSTRUCTIONS FOR USE:

### 1st step

Select and record the color of the patient's teeth using a color

scale and/or photograph before starting whitening.

### 2nd step

Perform the relative isolation with DSP dam (light-cured gingival barrier) covering the marginal gingiva and papillae with a layer of 3 to 5 mm wide and a maximum of 1 mm thick. The gingival barrier should cover approximately 0.5 to 1 mm of the dental surface.

### Notes:


Use a clinical mirror looking from incisal to cervical and observe for uncovered gingival tissue. If so, perform the correction. This step is crucial to avoid contact of Hydrogen Peroxide with the gum. Use a lip retractor to facilitate the application of the gingival barrier and whitener.

### 3rd step

Polymerize the light-cured gingival barrier DSP dam using 20 to 30 seconds of light curing for each group of 3 to 4 teeth. The gingival barrier that forms is rigid and insoluble, preventing possible irritation by aggressive products.

### 4th step

Remove the protective cap from the Phase A syringe, then do the same with the Phase B syringe. Connect the two syringes, screw until you feel their final locking. Inject the contents of the Phase A syringe into the Phase B syringe and repeat the same operation in the reverse direction until you have completed 20 times, that is, 10 movements in each direction smoothly. After mixing is complete, return the gel to the Phase A syringe. Carefully disconnect the syringes and adapt the applicator tip to the Phase A syringe. Dispense the gel onto the dental surfaces in sweeping movements. If air bubbles are noticed during application, remove them with the aid of a microapplicator or a spatula and eliminate the air bubble on a surface before continuing application of the product. Ensure that all surfaces in question are completely covered by the gel;

### 5th step

With the aid of a brush or spatula, fully cover the vestibular surface of the teeth to be whitened, including the interproximal teeth and extend a little on the incisal and occlusal faces. The gel layer should be between 0.5 and 1mm thick.

### 6th Step

It is suggested to maintain the product in contact with the dental surfaces for a period of 40 to 45 consecutive minutes, and there is no need for three consecutive changes every 15 minutes. This fact is due to the properties of the product in keeping its pH controlled and stable for a longer period of time after mixing. It is not necessary to apply a light source as an additional method to accelerate the lightening process. The product will act effectively even in the absence of photo-activation.

### 7th Step

At the end of the recommended time, suck the gel over the teeth with a vacuum cannula and clean them with a gauze.

### 8th Step

Remove the gingival barrier by highlighting it forward.

### 9th Step

Apply the DSP Desensitizer 2% for 10 minutes and then polish the teeth with DSP Diamond Polishing paste.

Check the final appearance after treatment. Perform a comparison with the start and end photo.

### Notes:

- Monitor the patient throughout the whitening with respect to dental sensitivity and possible points of peroxide irritation. If there are any of these cases, interrupt the process and perform a verification and intervention (see comments in the section Precautions and Side Effects).
- If the obtained whitening does not meet the expectation and if the patient does not present sensitivity or other contraindication, the reapplication of the product can be done for up to two more sessions. If more than one session is required, a minimum interval of 7 days should be given between sessions.

## PRECAUTIONS:

### Before Whitening

- Promote a detailed medical history.
- The patient should have a perfectly clean tooth surface and no calculus.
- The gums should be healthy without inflammation.
- Prophylaxis with water and pumice stone, not using sodium bicarbonate.
- Perform removal of calculus at least 2 days before whitening.
- If the patient has spontaneous sensitivity in the teeth, treat it before whitening.
- If the teeth to be whitened present gingival retraction, isolate the area with the DSP dam.

- Ensure that the teeth to be whitened do not contain cavities or infiltrations in the restorations.
- In case of teeth with small fracture or incisal wear, protect the area with DSP Dam before application of the whitening gel to avoid sensitivity.
- All gingival tissue adjacent to the teeth must be protected with the DSP dam.

## Upon Whitening

- Avoid contact of the gel with the soft tissues. If contact occurs, apply NEUTRALIZER.
- If the patient has sensitivity, the whitening can be interrupted after careful analysis by the professional. After application of the whitening gel, use the DSP Desensitizer 2% NaF and let it act for 3-5 minutes and polish with felt discs on the DENSITIZER.

## PRECAUTIONS AFTER WHITENING

- Do not smoke for at least 2 days.
- Do not ingest acid foods and soft drinks.
- For 24 hours do not ingest products containing dyes.
- Make use of white toothpaste and without sodium bicarbonate.

## CONTRAINDICATIONS:

- The product is not indicated for whitening in pregnant women, infants and patients under 15 years of age.
- Patients with periodontal pathologies.
- Patients with incomplete root formation.
- Patients with allergic reaction to the product components.

## GENERAL SAFETY INSTRUCTIONS:

- The product must be handled exclusively in the dental clinic.
- Wear gloves during the procedure.
- Wear safety goggles on the patient and the professional.
- Protect the lips with insulators such as petroleum jelly or similar.
- The contents of the peroxide bottles may be under pressure and/or may have leaked due to adverse storage conditions. Before opening the packaging, protect your hands and eyes and ensure that there is no leaked liquid inside.
- Avoid contact of the product with skin, mucous membranes or clothing.
- Do not ingest the product.
- In case of accidental contact with the skin or intraoral mucosa, the region should be cleaned with plenty of water. In these cases, the action of the whitening gel can be neutralized by covering the affected area with a few drops of the NEUTRALIZING solution.
- In case of accidental contact with eyes, rinse immediately with plenty of water and seek medical advice. Leftovers or residues from used materials should not be reused. Discard properly.

## SIDE EFFECTS:

Some patients may experience dental sensitivity during or after whitening. If the patient has sensitivity during treatment, discontinue whitening and treat the teeth with desensitizer DSP Desensitizer 2%. If the sensitivity persists and it is not possible to continue the whitening, schedule a new session with a minimum interval of 7 days. Assess whether there are any anomalies in the teeth that can be corrected to eliminate such sensitivity. There are cases where sensitivity appears only after treatment (1 or 2 hours after). In these cases, if the sensitivity is high, it is recommended to prescribe an analgesic and, eventually, an anti-inflammatory, according to the need of each case.

Whitening may enhance pre-existing areas of hypocalcification (hypoplasias) due to the different (more intense) whitening behavior of these areas. In some cases, these spots disappear due to the dehydration of the enamel after the whitening is complete. Fluoride treatment of these regions can also help in their disappearance. Discontinue use of the product in case of allergic reactions.

## CONSERVATION AND STORAGE:

- Keep the product in its packaging always closed.
- Protect from the incidence of direct sunlight.
- Store the product at temperatures from 0°C to 30°C.
- If they find malfunctions, request the replacement of the material.

## DATE OF MANUFACTURE: PACKING

EXPIRATION DATE: 2 YEARS AFTER THE DATE OF MANUFACTURE.

## WARNING:

Respect the expiration date. Never use the product if it has expired. Keep out of reach of children. For dental use only. Before using observe the appearance of the product, if you notice any anomalies, contact the manufacturer or a professional immediately. The manufacturer is not responsible for damages caused by other uses or by incorrect handling. This material has been manufactured for dental use only and must be handled according to the instructions for use. It is the responsibility of the user to prove, before the use of this product, that this material is compatible and suitable for the desired use, especially when this use is not indicated in these instructions for use. Data descriptions do not constitute any kind of guarantee and, therefore, do not have any link. Do not reuse empty packages. Do not

reuse empty packaging for other purposes. To dispose of the product, follow local legislation.

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**PROIBIDO REPROCESSAR**

**Prohibited Reprocess**  
Proibido Reprocessar  
Prohibido reprocessar

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