Force post

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EN | ENGLISH

PT | PORTUGUÊS

ES | ESPAÑOL

DSP
white

Force post

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ENGLISH

To access the IFU go to dspbiomedical.com/ifu with code: 0.9000

TRADE NAME: DSP ForcePost
TECHNICAL NAME: Fiberglass Pin

Please read all information of this instruction for use carefully before applying the product.

INSTRUCTIONS FOR USE PROFESSIONAL USE ONLY

Keep for consultation until the total consumption of the product.
Keep this instruction until there is no more interaction of the product with your patient.

CONTENT:

Forms of presentation:
Complete Kit (DSP ForcePost DC):
Contains 5 pins of each size (DC0.5, DC1, DC2, DC3 and DC2L) with their respective bits (25 pins and 5 bits).

Special Kit (DSP ForcePost DC and DCL):
Contains 5 pins of each size (DC0.5, DC0.5L, DC1, DCL and DC2L) with their respective bits (25 pins and 5 bits).

Introductory Kit (DSP ForcePost DC or DCL):
Contains 5 DSP ForcePost pins of a given size and their respective Drill bit.

Refill (DSP ForcePost DC or DCL):
Contains 05 ForcePost DSP pins of a given size.

Unit Refill (DSP ForcePost DC or DCL):
Contains 1 ForcePost DSP pins of a given size.

PRODUCT DESCRIPTION:

DSP ForcePost is a pin made of fiberglass composite and high mechanical strength epoxy resin, which acts as an intraradicular reinforcement of the dental structure. The DSP ForcePost has a silane load pretreatment and promotes retention for the definitive restorative material or filling core, in the case of indirect restorations.
The fiberglass posts DSP ForcePost DCL aim to meet cases where there is greater wear of the root canal and consequently greater mechanical effort in the cervical third of the anterior teeth.

CÓDIGO	DESCRIÇÃO	A	B	C
DSP.9142-1	Forcepost 0.5 DC	1.4	20	0.65
DSP.9181-1	Forcepost 0.5 DCL	1.8	18	0.65
DSP.9162-1	Forcepost 1 DC	1.6	20	0.85
DSP.9201-1	Forcepost 1 DCL	2.0	18	0.85
DSP.9182-1	Forcepost 2 DC	1.8	20	1.05
DSP.9221-1	Forcepost 2 DCL	2.2	18	1.05
DSP.9202-1	Forcepost 3 DC	2.0	20	1.25
DSP.9222-1	Forcepost 4 DC	2.2	20	1.45

PRODUCT INDICATION:

The indication of an intraradicular fiberglass pin is directly related to the need to restore an endodontically treated tooth to assist in the retention of the restorative material and distribute stresses imposed on the tooth. In the case of a tooth with great coronal destruction due to caries, trauma or extensive restoration, it is recommended to place an intraradicular pin to reinforce the coronal portion, thus minimizing the probability of fracture. It can be used with dual or light-cured adhesive systems due to its great capacity to transmit light.

COMPOSITION:

Pins (% in weight): Fiberglass: 80% - Epoxy Resin: 20%

Property	Specified
Elastic Modulus * (GPa)	30 to 50
Bending Strength * (MPa)	Minimum 1100

INSTRUCTIONS FOR USE:

1st step
Perform the radiographic examination (periapical radiography) and check the quality of the endodontics. Select the pin according to the diameter of the conduit, positioning the DSP ForcePost Selection Ruler or the pin on the tooth radiography;

2nd step
Isolate the operative field. Prioritize the use of the rubber dam;

3rd step
Remove the gutta percha from the channel using a drill or suitable instrument keeping 3 to 5 mm of gutta percha in the apical region, ensuring the sealing of the endodontic filling;

4th step
Prepare the root canal with the specific drill bits for each pin (number) for optimal adaptation;

5th step
Place the selected pin in position and x-ray to evaluate its adaptation and its positioning. The pin must be perfectly adapted to the preparation performed, not presenting spaces between the end of the pin and the gutta percha;

6th Step
Cut out the coronal excess of the pin if necessary and take it back in position to check the adjustment;

NOTE 1: If the conduit is too wide, opt for a DCL ForcePost DSP pin.

NOTE 2: Clean the DSP ForcePost pin with 70% alcohol or DSP Asepsis 2% Solution (DSP).
It is recommended to sterilize the pin in an autoclave (121°C and 15 psi) before cementing. The pin should be very dry when cementing.

Use a proper packaging during sterilization and observe the sterilization conditions specified. Conditions other than those specified may damage the product.


7th Step
Apply phosphoric acid gel (DSP Etch 37%) for 15 seconds inside the root canal and remaining tooth structure. Wash thoroughly and remove excess water with absorbent paper cones;

8th Step
Use the Microapplicator for the application of the adhesive system according to the manufacturer's instructions throughout the length of the conditioned root canal and light cure for 40 seconds;

9th Step
Apply silane (DSP Silane) on the surface of the entire pin and wait for it to dry;

Step 10
With a slow drill, take the resin cement (DSP Force Dual Cure) to the channel. If you prefer, smear the pin with cement and insert it into the conduit;

Step 11
Carefully insert the pin into the channel and proceed to light cure the cement by the surface and through the pin.

PRECAUTIONS AND CONTRAINDICATIONS:

The pins may have exposed tips of the glass fibers that compose them and these can pierce the skin. The cutting

of the pins to adjust the length can also expose perforating fibers. Wear goggles during the cutting of the pins and avoid cutting near the patient's face to prevent any possible accident or discomfort with cut residues. Avoid vacuuming the dust from the cutting of the pins. To handle the product always use latex gloves. The pins must be free of oily dirt, organic materials or any other contaminants at the time of cementation to avoid adhesion failure.

CONSERVATION AND STORAGE:

- Keep the product in its original packaging.
- Protect from the incidence of direct sunlight.
- Store the product at temperatures from 0°C to 30°C.
- If they find malfunctions, request the replacement of the material.

DATE OF MANUFACTURE: PACKING

EXPIRATION DATE: 5 YEARS AFTER THE DATE OF MANUFACTURE.

WARNING:

Respect the expiration date. Never use the product if it has expired.
Keep out of reach of children.
For dental use only. Before using, observe the appearance of the product, if you notice any anomalies, contact the manufacturer immediately. The manufacturer is not responsible for damage caused by improper use or improper handling.
This product has been manufactured for dental use only and must be handled according to the instructions for use. It is the professional's responsibility to prove, before using the product, whether this product is compatible and suitable for the desired use.
Do not reuse empty packaging for other purposes. To dispose of the product, follow local legislation.

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