



DSP PORCELAIN ETCH

MANUFACTURED AND DISTRIBUTED BY:



DSP Industrial LTDA | CNPJ: 03.960.018/0001-23

Rua Marechal Floriano Peixoto, 303 Ouro Verde II - Campo Largo - PR 83.606-290
- SAC 0800.600.88.66

Resp. Técnico: Célio G. Netzel
CREA-PR: 25412/D

Anvisa Registration: 80116980024

FOR PURCHASES, QUESTIONS, COMPLIMENTS, AND COMPLAINTS: 0800 600 8866

INSTRUCTIONS FOR USE | english

TRADE NAME: DSP Porcelain Etch
TECHNICAL NAME: Conditioning Acid
CÓDIGO: 0.1200

Please read all information of this instruction for use carefully before applying the product.

INSTRUCTIONS FOR USE PROFESSIONAL USE ONLY

Keep for consultation until the total consumption of the product. Keep this instruction until there is no more interaction of the product with your patient.

CONTENT:

- 1 acid conditioner syringes containing 3g with 1 applicator tip,
- 1 instruction manual for the professional and patient.

PRODUCT DESCRIPTION:

DSP Porcelain Etch is a hydrofluoric acid-based gel in concentrations of 5 or 10%, for the acid conditioning of porcelains, with the objective of increasing the adhesion of resinous materials.

PRODUCT INDICATION:

Surface conditioning of dental porcelain for adhesion in the process of restorations, such as feldspathic porcelain, ceramic glass and fixations of orthodontic brackets only on porcelain teeth.

COMPOSITION:

Hydrofluoric Acid 5 or 10%, Dye, Thickeners and Deionized Water.

INSTRUCTIONS FOR USE:

1st Feldspathic Porcelain Condition for 1-2 minutes with DSP Porcelain Etch 10% HF;

DSP PORCELAIN ETCH

2nd Glass ceramic reinforced with Leucite Condition for 1 minute with DSP Porcelain Etch 10% HF;

3rd Lithium Disilicate based glass ceramic Condition for 20-30 seconds with DSP Porcelain Etch 5% HF;

4rd Zirconia Non-conditionable.

PRECAUTIONS AND CONTRAINDICATIONS:

Hydrofluoric acid is corrosive. Avoid contact with enamel, dentin, eyes, skin and mucosa. It is recommended to use eye protection for the patient, dentist and assistant during its use. When applied to the patient's mouth, the oral tissues must be protected with Vaseline and, whenever possible, perform absolute isolation. In case of accidental contact with eyes, skin or mucosa, rinse thoroughly with water. If necessary, seek medical attention.

SIDE EFFECTS:

Contact with soft tissues and ingestion of the product causes irritation and fluoride poisoning. Discontinue use of the product in case of allergic reactions.

CONSERVATION AND STORAGE:

- Keep the product in its packaging always closed.
- Protect from the incidence of direct sunlight.
- Store the product at temperatures from 0°C to 30°C.
- If they find malfunctions, request the replacement of the material.














DATE OF MANUFACTURE: PACKING

EXPIRATION DATE: 2 YEARS AFTER THE DATE OF MANUFACTURE

WARNING:

Respect the expiration date. Never use the product if it has expired. Keep out of reach of children. For dental use only. Before using observe the appearance of the product, if you notice any anomalies, contact the manufacturer or a professional immediately. The manufacturer is not responsible for damages caused by other uses or

by incorrect handling. This material has been manufactured for dental use only and must be handled according to the instructions for use. It is the responsibility of the user to prove, before the use of this product, that this material is compatible and suitable for the desired use, especially when this use is not indicated in these instructions for use. Data descriptions do not constitute any kind of guarantee and, therefore, do not have any link. Do not reuse empty packages. Do not reuse empty packaging for other purposes. To dispose of the product, follow local legislation.

	Code of Product
	Number of the lot
	Manufacturer
	Date of manufacturing
	Representative in the community european
	Prohibited Reprocess
	Date of Durability Minimum
	Do not reuse
	Do not use if the package it's damaged
	Keep Dry
	Keep away from light
	Notification demanded by FDA for Business In the USA
	Reference to information join or attached