PORTUGUÊS

Nome comercial
MINI IMPLANTE ORTHOFIT DSP

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ORTHODONTIC IMPLANT

Presentation

The MINI IMPLANT ORTHOFIT DSP is an implantable medical device designed to be attached to the bone, assisting orthodontic movements. Are a single-use product, sterilized by gamma radiation and, after orthodontic treatment, it should be removed and immediately disposed must comply with the health legislation recommended by the Health Surveillance Agency for contaminated products. Reprocessing is forbidden for this product. It's commercialized sterilized and packed in a blister containing an unit of the product.

Technical Product Features

MINI IMPLANTE ORTHOFIT DSP is an implantable device, made with Titanium ASTM F136. Your body has a conical shape self-drilling, available in the di-

ameters of 1.4; 1.5 e 1.6 and length with variations of 7,00 a 11.5 mm. It has height variation on belt (neck brace) of 0.5; 2.0; 2.4 e 4.0 mm. Its head have features of mechanical retention with chambers and a perforation for orthodontic wire fixation.

Manufacturing Material

MINI IMPLANT ORTHOFIT DSP made of Titanium ASTM F136.

Indication

The ORTHOFIT DSP Minimplants are indicated to assist orthodontic movements such as:

- Mesialization of molar teeth;
- Intrusion of lower and upper incisors;
- Intrusion of molars;
- Retraction of anterior and inferior teeth;
- Verticalization of molars;
- Distortion of premolars;
- Retreatment of canines.

Instructions for use

- Select the length and diameter of the ORTHOFIT DSP mini implant according to the interdental space, as well as the desired orthodontic movement. The low, medium or long straps depend on the area of implantation. If there is a bone dip or defect, a long strap has more indication because it will prevent attached accessories (eg springs or elastic chains) from leaning against the mucosa.
- With the help of the Contra-angle Connection wrench, capture the ORTHOFIT DSP Mini Implant.
- Take to the receiving area and install the ORTHOFIT DSP mini implant as planned.
- Pre-drilling with drills is not required. The ORTHOFIT DSP mini implant is self-drilling. Maximum torque for setting - 10 Ncm. The estimated time of permanence of the mini implant after installation may vary from 6 months to 2 years, or according to the results obtained during the treatment.

Contraindications

- Generally, the contraindications associated with:
- No bone support for anchorage.
- Anatomical defects in bone structures.
- Anatomical defects in cavities that harbor the maxillary sinus and alveolar nerve.
- Bone quality and quantity.

Risk prevention, precautions, restrictions

- Single-use product - REPROCESSING IS FORBIDDEN;
- Do not use the product if the package is damaged or expired;
- Products should be used with sterile and clean surgical instruments;
- In the event of reuse, there are risks such as: infection, patient contamination with other chemicals.
- The misuse, abuse or excessive force can cause stripping, breakage or irreversible damage to the product;
- The use of these products requires a specific professional expertise, so DSP is not responsible for technical errors of incorrect product application or improper product handling by the professional.
- Attach the extra label supplied with the product to the patient's registration form to enable product traceability.

Product Disposal

After the residence time (which varies according to each case), the product must be removed and immediately disposed of in accordance with current environmental regulations

Warning

The product is marketed sterile. It is not recommended to put it in an oven or autoclave.

Storage and transportation

Keep in a cool place away from moisture and sunlight

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